

For Research Use Only

TBC1D10C Polyclonal antibody, PBS Only

Catalog Number: 34255-1-PBS



Basic Information

Catalog Number:

34255-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40608

GenBank Accession Number:

BC062999

GeneID (NCBI):

374403

UNIPROT ID:

Q8IV04

Full Name:

TBC1 domain family, member 10C

Calculated MW:

446 aa, 50 kDa

Observed MW:

47 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 34255-1-AP (TBC1D10C antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 34255-1-PBS in a different storage buffer formulation.