

For Research Use Only

# GFP Monoclonal antibody (3H9)



[www.ptglab.com](http://www.ptglab.com)

Catalog Number: 3h9-2

## Basic Information

**Catalog Number:**  
3h9-2

**Applications:**  
WB

**Host/IsoType:**  
Rat/IgG2a

**Conjugate:**  
Unconjugated

**Type:**  
Primary Antibody

**RRID:**  
AB\_10773374

## Description

Rat monoclonal antibody [3H9] to Green Fluorescent Proteins

## Specificity/Target

The antibody recognizes green fluorescent proteins derived from *Aequorea victoria* GFP, such as eGFP, wtGFP, YFP, or CFP.

## Recommended Dilution

Western Blot: 1:1,000

## Storage

**Storage:**

Shipped on ice. Store at -20°C/-4°F. Stable for 1 year after shipment.

**Storage Buffer:**


PBS, 50% glycerol, preservative: 0.09% sodium azide.

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
1(312) 455-8498 (outside USA) W: [www.ptglab.com](http://www.ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Dilution series of EGFP in HEK 293T cell lysate  
Primary antibody: [3H9] 1:1,000 Secondary  
antibody: goat anti-rat Alexa Fluor® 488 1:1,000