

For Research Use Only

Ccr5 Polyclonal antibody

Catalog Number: 50650-1-AP

1 Publications



Basic Information

Catalog Number: 50650-1-AP	GenBank Accession Number: BC103587	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 147 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 12774	
Source: Rabbit	Full Name: chemokine (C-C motif) receptor 5	
Isotype: IgG	Calculated MW: 40 kDa	
Immunogen Catalog Number: AG0577		

Applications

Tested Applications:
ELISA

Cited Applications:
IF, IHC

Species Specificity:
human, mouse, rat

Cited Species:
monkey

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Julian B Hattler	37366075	Brain Behav	IF,IHC

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data