

For Research Use Only

ICMT Polyclonal antibody

Catalog Number: 51001-2-AP

Featured Product

5 Publications



Basic Information

Catalog Number:

51001-2-AP

Size:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC028168

GeneID (NCBI):

23463

UNIPROT ID:

O60725

Full Name:

isoprenylcysteine carboxyl methyltransferase

Calculated MW:

284 aa, 32 kDa

Observed MW:

30-34 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Helen Court	29051265	J Cell Biol	WB
Carla Borini Etichetti	30655292	J Biol Chem	WB
Ian M Ahearn	33579760	Life Sci Alliance	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

