

For Research Use Only

ICMT Polyclonal antibody

Catalog Number: 51001-2-AP

Featured Product

5 Publications



Basic Information

Catalog Number: 51001-2-AP	GenBank Accession Number: BC028168	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 23463	
Source: Rabbit	Full Name: isoprenylcysteine carboxyl methyltransferase	
Isotype: IgG	Calculated MW: 284 aa, 32 kDa	
	Observed MW: 30-34 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Cited Species:
human, mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Helen Court	29051265	J Cell Biol	WB
Carla Borini Etichetti	30655292	J Biol Chem	WB
Ian M Ahearn	33579760	Life Sci Alliance	WB

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data