

For Research Use Only

Evi2a Polyclonal antibody

Catalog Number: 51014-2-AP



Basic Information

Catalog Number:

51014-2-AP

Size:

150ul , Concentration: 250 ug/ml by Nanodrop and 250 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0478

GenBank Accession Number:

BC038124

GeneID (NCBI):

14017

UNIPROT ID:

P20934

Full Name:

ecotropic viral integration site 2a

Calculated MW:

24 kDa

Observed MW:

26 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : mouse liver tissue,

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 51014-2-AP (Evi2a antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.