

For Research Use Only

Prnd Polyclonal antibody

Catalog Number: 51026-2-AP



Basic Information

Catalog Number:

51026-2-AP

Size:

150ul , Concentration: 133 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0731

GenBank Accession Number:

BC025140

GeneID (NCBI):

26434

UNIPROT ID:

Q9QUG3

Full Name:

prion protein dublet

Calculated MW:

21 kDa

Observed MW:

34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : human skeletal muscle tissue,

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 51026-2-AP (Prnd antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.