

For Research Use Only

Harmonin Polyclonal antibody

Catalog Number: 51028-2-AP

1 Publications



Basic Information

Catalog Number:

51028-2-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0513

GenBank Accession Number:

BC016057

GeneID (NCBI):

10083

UNIPROT ID:

Q9Y6N9

Full Name:

Usher syndrome 1C (autosomal recessive, severe)

Calculated MW:

552 aa, 62 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IHC

Species Specificity:

human, mouse

Cited Species:

mouse

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-----------------------------|-----------|------------|-------------|
| Grimsley-Myers Cynthia M CM | 20016102 | J Neurosci | IHC |

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

