

For Research Use Only

MAP4K4 Polyclonal antibody

Catalog Number: 51088-1-AP **1 Publications**



Basic Information

Catalog Number: 51088-1-AP	GenBank Accession Number: BC156887	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 287 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 9448	
Source: Rabbit	Full Name: mitogen-activated protein kinase kinase kinase 4	
Isotype: IgG	Calculated MW: 142 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IHC, WB

Species Specificity:
human

Cited Species:
rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Lifeng Ding	37647682	Int Immunopharmacol	WB,IHC

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data