

For Research Use Only

VSP15a Polyclonal antibody

Catalog Number: 51090-1-AP



Basic Information

Catalog Number:

51090-1-AP

GenBank Accession Number:

CAC21478

Purification Method:

Antigen affinity purification

Size:

150ul , Concentration: 900 ug/ml by Nanodrop and 133 ug/ml by Bradford method using BSA as the standard;

GeneID (NCBI):

2543600

Full Name:

SPAPB1A10.05

Source:

Rabbit

Calculated MW:

32 kDa

Isotype:

IgG

Immunogen Catalog Number:

AG0724

Applications

Tested Applications:

ELISA

Species Specificity:

yeast

Background Information

This antibody recognize the SPAPB1A10 protein of Yeast.

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

