

For Research Use Only

Cxcl12 Polyclonal antibody

Catalog Number: 51091-1-AP

3 Publications



Basic Information

Catalog Number:

51091-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0474

GenBank Accession Number:

BC020297

GeneID (NCBI):

20315

UNIPROT ID:

P40224

Full Name:

chemokine (C-X-C motif) ligand 12

Calculated MW:

11 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IHC, IF

Species Specificity:

human

Cited Species:

human, rabbit

Notable Publications

Author	Pubmed ID	Journal	Application
Chunyan Xiang	38662111	Clin Exp Med	IHC
Toru Ichiseki	37443763	Cells	IF
Raimundo Fernandes de Araújo Júnior	36857852	Transl Oncol	IHC

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

