

For Research Use Only

CYP302a1 (mosquito) Polyclonal antibody



Catalog Number: 51097-1-AP

Basic Information

Catalog Number: 51097-1-AP	GenBank Accession Number: AY947549	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 160 µg/ml by Nanodrop and 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): Full Name: Calculated MW: 59 kDa	
Source: Rabbit		
Isotype: IgG		
Immunogen Catalog Number: AG0761		

Applications

Tested Applications:
ELISA

Species Specificity:
mosquito

Background Information

CYP302a1, also named as disembodied (dib-Bm), belongs to the cytochrome P450 family. This antibody detects the CYP302a1 protein of mosquito.

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data