

For Research Use Only

Collagen Type II Polyclonal antibody



Catalog Number: 51104-1-AP

9 Publications

Basic Information

Catalog Number: 51104-1-AP	GenBank Accession Number: NM_001844	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 147 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 1280	
Source: Rabbit	Full Name: collagen, type II, alpha 1	
Isotype: IgG	Calculated MW: 142 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IF, IHC, WB

Species Specificity:
human

Cited Species:
human, rat, mouse, rabbit

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhao Zhang	36125665	Apoptosis	WB,IHC
Yue Li	36459925	Int Immunopharmacol	WB
Zhilei Hu	36105191	Front Pharmacol	WB

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data