

For Research Use Only

Collagen Type II Polyclonal antibody

Catalog Number: 51104-1-AP

19 Publications



Basic Information

Catalog Number:

51104-1-AP

Size:

150ul, Concentration: 147 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001844

GeneID (NCBI):

1280

UNIPROT ID:

P02458

Full Name:

collagen, type II, alpha 1

Calculated MW:

142 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human

Cited Species:

human, mouse, rat, pig, rabbit

Notable Publications

Author	Pubmed ID	Journal	Application
Zhao Zhang	36125665	Apoptosis	WB, IHC
Yue Li	36459925	Int Immunopharmacol	WB
Zhilei Hu	36105191	Front Pharmacol	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

