

For Research Use Only

# Ccl8 Polyclonal antibody

Catalog Number: 51110-1-AP



## Basic Information

Catalog Number: 51110-1-AP	GenBank Accession Number: BC117101	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 20307	
Source: Rabbit	Full Name: chemokine (C-C motif) ligand 8	
Isotype: IgG	Calculated MW: 11 kDa	
Immunogen Catalog Number: AG0466		

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -20°C.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data