

For Research Use Only

Tbp Polyclonal antibody

Catalog Number: 51128-1-AP **1 Publications**



Basic Information

Catalog Number:

51128-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC016476

GeneID (NCBI):

21374

UNIPROT ID:

P29037

Full Name:

TATA box binding protein

Calculated MW:

28 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Notable Publications

Author	Pubmed ID	Journal	Application
Pan Sun	40919497	Stem Cells Int	

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data