

For Research Use Only

# Vpr Polyclonal antibody

Catalog Number: 51143-1-AP

Featured Product

10 Publications



## Basic Information

Catalog Number:

51143-1-AP

Size:

150ul, Concentration: 800 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0718

GenBank Accession Number:

NP\_057852

GeneID (NCBI):

Full Name:

Calculated MW:

11 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

immunodeficiency virus 1

Cited Species:

human, virus

## Background Information

This antibody recognize the VPR gene of Human immunodeficiency virus 1.

## Notable Publications

Author	Pubmed ID	Journal	Application
Qiong Zhang	34545078	Nat Commun	WB
Sang Yoon Park	25352594	J Biol Chem	WB
Nicolas Vabret	35789859	iScience	WB

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

