

For Research Use Only

Vpr Polyclonal antibody

Catalog Number: 51143-1-AP

Featured Product

12 Publications



Basic Information

Catalog Number:

51143-1-AP

Size:

150ul, Concentration: 800 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0718

GenBank Accession Number:

NP_057852

GeneID (NCBI):

Full Name:

Calculated MW:

11 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

immunodeficiency virus 1

Cited Species:

human, virus

Background Information

This antibody recognize the VPR gene of Human immunodeficiency virus 1.

Notable Publications

Author	Pubmed ID	Journal	Application
Qiong Zhang	34545078	Nat Commun	WB
Sang Yoon Park	25352594	J Biol Chem	WB
Nicolas Vabret	35789859	iScience	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

