

For Research Use Only

MTCO1 Polyclonal antibody

Catalog Number: 55071-1-AP

4 Publications



Basic Information

Catalog Number:

55071-1-AP

Size:

150ul, Concentration: 287 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

YP_003024028

GeneID (NCBI):

4512

Full Name:

cytochrome c oxidase I

Calculated MW:

57 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Jing Zhang	26492917	EMBO J	WB
Hao-Peng Lin	27165101	Biomed Microdevices	WB
Xiaoting Liu	35876654	EMBO Rep	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

