

For Research Use Only

# Integrin Alpha 7 Polyclonal antibody

Catalog Number: 55073-1-AP



---

## Basic Information

Catalog Number:	55073-1-AP	GenBank Accession Number:	NM_002206	Purification Method:	Antigen affinity purification
Size:	150ul	GeneID (NCBI):	3679		
Source:	Rabbit	Full Name:	integrin, alpha 7		
Isotype:	IgG	Calculated MW:	129 kDa		

---

## Applications

Tested Applications:  
ELISA  
Species Specificity:  
human

---

## Background Information

---

### Storage

Storage:  
Store at -20°C.  
Storage Buffer:  
0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

---

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data