

For Research Use Only

GRM5 Polyclonal antibody

Catalog Number: 55081-1-AP



Basic Information

Catalog Number:	55081-1-AP	GenBank Accession Number:	NM_000842	Purification Method:	Antigen affinity purification
Size:	150ul	GeneID (NCBI):	2915		
Source:	Rabbit	Full Name:	glutamate receptor, metabotropic 5		
Isotype:	IgG	Calculated MW:	132 kDa		

Applications

Tested Applications:
ELISA
Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C.
Storage Buffer:
0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data