

For Research Use Only

CSPG5 Polyclonal antibody

Catalog Number: 55086-1-AP



Basic Information

Catalog Number:

55086-1-AP

Size:

150ul, Concentration: 240 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006574

GeneID (NCBI):

10675

UNIPROT ID:

O95196

Full Name:

chondroitin sulfate proteoglycan 5 (neuroglycan C)

Calculated MW:

60 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

CSPG5, also named as CALEB and NGC, may function as a growth and differentiation factor involved in neuritogenesis. It may induce ERBB3 activation. The antibody is specific to CSPG5.

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

