

For Research Use Only

# GRIN2D Polyclonal antibody

Catalog Number: 55094-1-AP



## Basic Information

**Catalog Number:**

55094-1-AP

**Size:**

150ul , Concentration: 300 µg/ml by Bradford method using BSA as the standard;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_000836

**GeneID (NCBI):**

2906

**UNIPROT ID:**

O15399

**Full Name:**

glutamate receptor, ionotropic, N-methyl D-aspartate 2D

**Calculated MW:**

144 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human

## Storage

**Storage:**

Store at -20°C.

**Storage Buffer:**

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

