

For Research Use Only

GRIN2D Polyclonal antibody

Catalog Number: 55094-1-AP



Basic Information

Catalog Number: 55094-1-AP	GenBank Accession Number: NM_000836	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 2906	
Source: Rabbit	Full Name: glutamate receptor, ionotropic, N-methyl D-aspartate 2D	
Isotype: IgG	Calculated MW: 144 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C.

Storage Buffer:
0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data