

For Research Use Only

RCOR1-Specific Polyclonal antibody

Catalog Number: 55174-1-AP



Basic Information

Catalog Number:
55174-1-AP

Size:
150ul

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
NM_015156

GeneID (NCBI):
23186

Full Name:
REST corepressor 1

Calculated MW:
53 kDa

Purification Method:
Antigen affinity purification

Recommended Dilutions:
IP 0.5-4.0 ug for IP and 1:500-1:1000
for WB

Applications

Tested Applications:
IP, ELISA

Species Specificity:
human, mouse

Positive Controls:

IP: mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data