

For Research Use Only

# GOLGA4 Polyclonal antibody, PBS Only

Catalog Number: 55207-1-PBS



## Basic Information

**Catalog Number:**

55207-1-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_002078

**GeneID (NCBI):**

2803

**UNIPROT ID:**

Q13439

**Full Name:**

golgi autoantigen, golgin subfamily a, 4

**Calculated MW:**

161 kDa

**Observed MW:**

260-270 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

WB, Indirect ELISA

**Species Specificity:**

human, mouse

## Background Information

GOLGA4, also known as Golgi matrix protein, is a member of the Golgin family and is mainly localized in the Golgi apparatus. It is involved in vesicle transport at the Golgi level and may play a role in the transport of GPI-anchored proteins from the Golgi network. GOLGA4 is expressed in a variety of tissues, especially in the testis, but its specific function in spermatogenesis and male fertility is not fully understood. Studies have shown that GOLGA4 knockout male mice are still fertile, and the testis and sperm morphology are normal, indicating that GOLGA4 is not essential for mouse spermatogenesis and male fertility (PMID: 32736686).

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 55207-1-AP (GOLGA4 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 55207-1-PBS in a different storage buffer formulation.