

For Research Use Only

DOT1L Polyclonal antibody

Catalog Number: 55315-1-AP



Basic Information

Catalog Number: 55315-1-AP	GenBank Accession Number: NM_032482	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 353 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 84444	
Source: Rabbit	Full Name: DOT1-like, histone H3 methyltransferase (S. cerevisiae)	
Isotype: IgG	Calculated MW: 185 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data