

For Research Use Only

PPP1R3F Polyclonal antibody

Catalog Number: 55386-1-AP



Basic Information

Catalog Number: 55386-1-AP	GenBank Accession Number: NM_033215	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop and 760 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 89801	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: protein phosphatase 1, regulatory (inhibitor) subunit 3F	
Isotype: IgG	Calculated MW: 83 kDa	
	Observed MW: 83 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : SMMC-7721 cells,
Species Specificity: human, mouse	

Background Information

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 55386-1-AP (PPP1R3F Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.