

For Research Use Only

# TGF Beta Receptor I Polyclonal antibody

Catalog Number: 55391-1-AP

4 Publications

## Basic Information

### Catalog Number:

55391-1-AP

### Size:

150ul , Concentration: 320 µg/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### GenBank Accession Number:

NM\_004612

### GeneID (NCBI):

7046

### UNIPROT ID:

P36897

### Full Name:

transforming growth factor, beta receptor 1

### Calculated MW:

56 kDa

### Observed MW:

56 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

ELISA

### Cited Applications:

WB, IHC

### Species Specificity:

human, mouse

### Cited Species:

human, mouse

## Notable Publications

Author	Pubmed ID	Journal	Application
Na Tang	28942050	Exp Parasitol	WB
Jingjing Wu	29077259	J Cell Mol Med	WB
Wanming He	30742067	Oncogene	WB

## Storage

### Storage:

Store at -20°C.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

