

For Research Use Only

BCL11B Polyclonal antibody, PBS Only

Catalog Number: 55414-1-PBS



Basic Information

Catalog Number:

55414-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_022898

GeneID (NCBI):

64919

UNIPROT ID:

Q9C0K0

Full Name:

B-cell CLL/lymphoma 11B (zinc finger protein)

Calculated MW:

96 kDa

Observed MW:

116 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, IP, Indirect ELISA

Species Specificity:

human, mouse

Background Information

BCL11B, also named CTIP2, RIT1, CTIP-2, ZNF856B, and hRIT1-alpha, is a tumor-suppressor protein involved in T-cell lymphomas. It may function on the P53-signaling pathway. BCL11B is a key regulator of both differentiation and survival during thymocyte development. It depresses transcription through direct, TF耦UP2-independent binding to a GC-rich response element. This antibody is specific to BCL11B. The calculated molecular weight of BCL11B is 95 kDa but the modified BCL11B protein is about 116 kDa.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 55414-1-AP (BCL11B antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 55414-1-PBS in a different storage buffer formulation.



IP result of anti-BCL11B (IP:55414-1-AP, 4ug; Detection:55414-1-AP 1:5000) with Jurkat cells lysate 1320 ug. This data was developed using the same antibody clone with 55414-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 55414-1-AP (BCL11B antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 55414-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using 55414-1-AP (BCL11B antibody) at dilution of 1:50 and CoraLite488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 55414-1-PBS in a different storage buffer formulation.