

For Research Use Only

RFP Monoclonal antibody (5F8)



www.ptglab.com

Catalog Number: 5f8-2

2 Publications

Basic Information

Catalog Number:
5f8-2

Clone No.:
5F8

Applications:
IF, ELISA

Host/IsoType:
Rat/IgG2a

Conjugate:
Unconjugated

Type:
Primary Antibody

RRID:
AB_2336064

Description

Rat monoclonal antibody [5F8] to Red Fluorescent Proteins

Specificity/Target

Tested on DsRed, mRFP, mCherry, mPlum, mRFPruby, mScarlet, tdTomato, mKate2, mOrange.

Recommended Dilution

IF: recommended starting dilution 1:1,000
A titration from 1:200 up to 1:2,000 is recommended.

Storage

Storage:
Shipped on ice. Store at -20°C/-4°F. Stable for 1 year after shipment.

Storage Buffer:
PBS, 50% glycerol, preservative: 0.09% sodium azide.

For technical support and original validation data for this product please contact

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data