

For Research Use Only

Cdc20 Monoclonal antibody

Catalog Number: 60010-1-Ig



Basic Information

Catalog Number:

60010-1-Ig

Size:

150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG2c

Immunogen Catalog Number:

AG0295

GenBank Accession Number:

BC001088

GeneID (NCBI):

991

UNIPROT ID:

Q12834

Full Name:

cell division cycle 20 homolog (S. cerevisiae)

Calculated MW:

55 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

1E8A8

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

