

For Research Use Only

# GATA1 Monoclonal antibody

Catalog Number: 60011-2-Ig



## Basic Information

<b>Catalog Number:</b> 60011-2-Ig	<b>GenBank Accession Number:</b> BC009797	<b>Purification Method:</b> Caprylic acid/ammonium sulfate precipitation
<b>Size:</b> 150ul	<b>GeneID (NCBI):</b> 2623	<b>CloneNo.:</b> 5B1A7
<b>Source:</b> Mouse	<b>Full Name:</b> GATA binding protein 1 (globin transcription factor 1)	
<b>Isotype:</b> IgG1	<b>Calculated MW:</b> 43 kDa	
<b>Immunogen Catalog Number:</b> AG1350		

## Applications

**Tested Applications:**  
WB, ELISA

**Species Specificity:**  
mouse

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com      W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data