

For Research Use Only

IFNAR2 Monoclonal antibody

Catalog Number: 60013-1-Ig



Basic Information

Catalog Number: 60013-1-Ig	GenBank Accession Number: BC002793	Purification Method: Caprylic acid/ammonium sulfate precipitation
Size: 150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 3455	CloneNo.: 1D4D9
Source: Mouse	Full Name: interferon (alpha, beta and omega) receptor 2	
Isotype: IgA	Calculated MW: 331 aa, 37 kDa	
Immunogen Catalog Number: AG0191	Observed MW: 58-150 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data