

For Research Use Only

# IFNAR2 Monoclonal antibody

Catalog Number: 60013-1-Ig



## Basic Information

Catalog Number:

60013-1-Ig

Size:

150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgA

Immunogen Catalog Number:

AG0191

GenBank Accession Number:

BC002793

GeneID (NCBI):

3455

UNIPROT ID:

P48551

Full Name:

interferon (alpha, beta and omega) receptor 2

Calculated MW:

331 aa, 37 kDa

Observed MW:

58~150 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

1D4D9

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

