

For Research Use Only

IFNAR2 Monoclonal antibody

Catalog Number: 60013-2-Ig



Basic Information

Catalog Number:

60013-2-Ig

Size:

150ul

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG0191

GenBank Accession Number:

BC002793

GeneID (NCBI):

3455

UNIPROT ID:

P48551

Full Name:

interferon (alpha, beta and omega)
receptor 2

Calculated MW:

331 aa, 37 kDa

Purification Method:

Caprylic acid/ammonium sulfate
precipitation

CloneNo.:

4G8G1

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

