

For Research Use Only

TIRAP Monoclonal antibody

Catalog Number: 60016-1-Ig



Basic Information

Catalog Number: 60016-1-Ig	GenBank Accession Number: BC032474	Purification Method: Protein G purification
Size: 150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 114609	CloneNo.: 4H5E12
Source: Mouse	UNIPROT ID: P58753	
Isotype: IgG1	Full Name: toll-interleukin 1 receptor (TIR) domain containing adaptor protein	
Immunogen Catalog Number: AG0771	Calculated MW: 256 aa, 28 kDa	
	Observed MW: 35 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

