

For Research Use Only

# RARRES2 Monoclonal antibody

Catalog Number: 60018-1-Ig **1 Publications**



## Basic Information

Catalog Number:	60018-1-Ig	GenBank Accession Number:	BC000069	Purification Method:	Caprylic acid/ammonium sulfate precipitation
Size:	150ul	GeneID (NCBI):	5919	CloneNo.:	2C4D9
Source:	Mouse	Full Name:	retinoic acid receptor responder (tazarotene induced) 2		
Isotype:	IgG1	Calculated MW:	163 aa, 19 kDa		
Immunogen Catalog Number:	AG0261				

## Applications

Tested Applications:  
ELISA

Cited Applications:  
WB

Species Specificity:  
human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Xiao-Hui Song	34957549	J Cell Physiol	WB

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data