

For Research Use Only

# RARRES2 Monoclonal antibody

Catalog Number: 60018-1-Ig **2 Publications**



## Basic Information

Catalog Number:

60018-1-Ig

Size:

150ul

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0261

GenBank Accession Number:

BC000069

GeneID (NCBI):

5919

UNIPROT ID:

Q99969

Full Name:

retinoic acid receptor responder  
(tazarotene induced) 2

Calculated MW:

163 aa, 19 kDa

Purification Method:

Caprylic acid/ammonium sulfate  
precipitation

CloneNo.:

2C4D9

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Xiao-Hui Song	34957549	J Cell Physiol	WB
Hideaki Nitta	38999437	J Clin Med	IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

