

For Research Use Only

# PRIM1 Monoclonal antibody

Catalog Number: 60021-1-Ig



## Basic Information

<b>Catalog Number:</b> 60021-1-Ig	<b>GenBank Accession Number:</b> BC005266	<b>Purification Method:</b> Caprylic acid/ammonium sulfate precipitation
<b>Size:</b> 150ul	<b>GeneID (NCBI):</b> 5557	<b>CloneNo.:</b> 5H2C8
<b>Source:</b> Mouse	<b>Full Name:</b> primase, DNA, polypeptide 1 (49kDa)	
<b>Isotype:</b> IgG2b	<b>Calculated MW:</b> 50 kDa	
<b>Immunogen Catalog Number:</b> AG1124		

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com      W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data