

For Research Use Only

# S100A11 Monoclonal antibody, PBS Only



Catalog Number: 60024-1-PBS

Featured Product

## Basic Information

Catalog Number:

60024-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0343

GenBank Accession Number:

BC001410

GeneID (NCBI):

6282

UNIPROT ID:

P31949

Full Name:

S100 calcium binding protein A11

Calculated MW:

105 aa, 12 kDa

Observed MW:

12 kDa

Purification Method:

Protein G purification

CloneNo.:

1A3F1

## Applications

Tested Applications:

WB, IF, IHC, Indirect ELISA

Species Specificity:

human

## Background Information

S100A11, also named as MLN 70, Calgizzarin and S100C, belongs to the S-100 family. It facilitates the differentiation and the cornification of keratinocytes. The naming of S100A11 systematically arises from its membership of the S100 protein family, named after their solubility in 100% saturated ammonium sulfate solution. First discovered in 1989, S100A11 has since been proposed to have direct biological functions in an assortment of physiological processes such as endo- and exocytosis, regulation of enzyme activity, cell growth and regulation, apoptosis and inflammation.

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




T-47D cells were subjected to SDS PAGE followed by western blot with 60024-1-Ig (S100A11 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60024-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 60024-1-Ig (S100A11 antibody) at dilution of 1:1000 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 60024-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 60024-1-Ig (S100A11 antibody) at dilution of 1:1000 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 60024-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed BxPC-3 cells using S100A11 antibody (60024-1-Ig, Clone: 1A3F1 ) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 60024-1-PBS in a different storage buffer formulation.