

For Research Use Only

ATF4 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60035-3-PBS



Basic Information

Catalog Number:

60035-3-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1279

GenBank Accession Number:

BC022088

GeneID (NCBI):

468

UNIPROT ID:

P18848

Full Name:

activating transcription factor 4 (tax-
responsive enhancer element B67)

Calculated MW:

39 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E12C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

60035-3-PBS targets ATF4 as part of a matched antibody pair:

MP51019-2: 60035-2-PBS capture and 60035-3-PBS detection (validated in Cytometric bead array)

MP51019-4: 60035-3-PBS capture and 60706-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51019-2, ATF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60035-2-PBS. Detection antibody: 60035-3-PBS. Standard:Ag16665. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51019-4, ATF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60035-3-PBS. Detection antibody: 60706-2-PBS. Standard:Ag16665. Range: 1.563-100 ng/mL