

For Research Use Only

MAGEC2 Monoclonal antibody

Catalog Number: 60040-1-Ig



Basic Information

Catalog Number:

60040-1-Ig

Size:

150ul, Concentration: 1000 µg/ml by Bradford method using BSA as the standard;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG0418

GenBank Accession Number:

BC005891

GeneID (NCBI):

51438

UNIPROT ID:

Q9UBF1

Full Name:

melanoma antigen family C, 2

Calculated MW:

373 aa, 41 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

2H5A12

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Melanoma antigen family C, 2(MAGEC2, synonyms: CT10, HCA587, MAGEE1, MGC13377) is related to members of the MAGEC gene family. MAGEC2 is not expressed in normal tissues, except for testis, and is expressed in tumors of various histological types. This gene and other MAGEC genes are clustered on chromosome Xq26-q27.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

