

For Research Use Only

OPN3 Monoclonal antibody

Catalog Number: 60090-1-Ig



Basic Information

Catalog Number: 60090-1-Ig	GenBank Accession Number: BC036773	Purification Method: Caprylic Acid/Ammonium Sulfate Precipitation
Size: 150ul , Concentration: 500 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 23596	CloneNo.: 1E10D10
Source: Mouse	Full Name: opsin 3	
Isotype: IgG1	Calculated MW: 402 aa, 45 kDa	
Immunogen Catalog Number: AG5075	Observed MW: 45 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data