

For Research Use Only

PFN2 Monoclonal antibody

Catalog Number: 60094-1-Ig **3 Publications**



Basic Information

Catalog Number: 60094-1-Ig	GenBank Accession Number: BC018049	Purification Method: Protein A purification
Size: 150ul	GeneID (NCBI): 5217	CloneNo.: 1D9C4
Source: Mouse	UNIPROT ID: P35080	
Isotype: IgM	Full Name: profilin 2	
Immunogen Catalog Number: AG5043	Calculated MW: 140 aa, 15 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB, IF

Species Specificity:
human, mouse

Cited Species:
human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Min Yu	29773903	Oncogene	WB
Sonia Afroz	27136678	Elife	IF
Xiaobin Cui	29238898	Mol Cell Biochem	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

