

For Research Use Only

APOBEC3G Monoclonal antibody

Catalog Number: 60100-2-Ig



Basic Information

Catalog Number: 60100-2-Ig	GenBank Accession Number: BC024268	Purification Method: Protein A purification
Size: 150ul	GeneID (NCBI): 60489	CloneNo.: 5G10A8
Source: Mouse	Full Name: apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3G	
Isotype: IgG2b	Calculated MW: 46 kDa	
Immunogen Catalog Number: AG7770		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

APOBEC3G (apolipoprotein B mRNA-editing enzyme, catalytic polypeptide-like 3G, also known as CEM15), a member of the APOBEC family of cytidine deaminase-editing enzymes, is a potent DNA mutator that can be packaged into budding retroviruses. HA-tagged APOBEC3G can be expressed in 293T cells, which do not express APOBEC3G endogenously, in addition 22- and 24-kDa proteolytic degradation fragments can be detected by western blot (PMID:14672928). It has 2 isoforms produced by alternative splicing.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data