

For Research Use Only

# TWEAKR Monoclonal antibody

Catalog Number: 60104-1-Ig



## Basic Information

Catalog Number:

60104-1-Ig

Size:

150ul

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG2689

GenBank Accession Number:

BC002718

GeneID (NCBI):

51330

UNIPROT ID:

Q9NP84

Full Name:

tumor necrosis factor receptor  
superfamily, member 12A

Calculated MW:

129 aa, 14 kDa

Purification Method:

Caprylic Acid/Ammonium Sulfate  
Precipitation

CloneNo.:

2E10C10

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human heart tissue were subjected to SDS PAGE followed by western blot with 60104-1-Ig (TWEAKR antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.