

For Research Use Only

FGFR2 Monoclonal antibody

Catalog Number: 60106-1-Ig



Basic Information

Catalog Number: 60106-1-Ig	GenBank Accession Number: BC039243	Purification Method: Protein A purification
Size: 150ul , Concentration: 2400 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 2263	CloneNo.: 5E1D5
Source: Mouse	UNIPROT ID: P21802	
Isotype: IgG2b	Full Name: fibroblast growth factor receptor 2	
Immunogen Catalog Number: AG3701	Calculated MW: 709 aa, 79 kDa	
	Observed MW: 130 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

