

For Research Use Only

FXYD6 Monoclonal antibody

Catalog Number: 60111-1-Ig



Basic Information

Catalog Number:

60111-1-Ig

Size:

150ul

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG8538

GenBank Accession Number:

BC018652

GeneID (NCBI):

53826

UNIPROT ID:

Q9HQ3

Full Name:

FXD domain containing ion transport
regulator 6

Calculated MW:

95 aa, 11 kDa

Purification Method:

Caprylic Acid/Ammonium Sulfate
Precipitation

CloneNo.:

6H10F4

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

