

For Research Use Only

# FXVD6 Monoclonal antibody

Catalog Number: 60111-1-Ig



## Basic Information

Catalog Number: 60111-1-Ig	GenBank Accession Number: BC018652	Purification Method: Caprylic Acid/Ammonium Sulfate Precipitation
Size: 150ul	GeneID (NCBI): 53826	CloneNo.: 6H10F4
Source: Mouse	Full Name: FXVD domain containing ion transport regulator 6	
Isotype: IgG2a	Calculated MW: 95 aa, 11 kDa	
Immunogen Catalog Number: AG8538		

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data